



2701 Bartram Road • Bristol, PA 19007
800-840-9041 • 215-788-2700
www.gentell.com • info@gentell.com

Comfortell® Dressing

4" x 4" (sterile) (2.5"x2.5" pad with 3/4" adhesive border) 75 dressings / box (75 dressings) GEN-12400 NDC# 61554-124-04	6" x 6" (sterile) (4.5"x4.5" pad with 3/4" adhesive border) 10/carton 5cartons / case (50 dressings) GEN-12600 NDC# 61554-126-06
-----------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------

- Sterile
- Water resistant border enables normal bathing
- Four distinct layers absorb mild to moderate exudate



Gentell Comfortell is a composite wound dressing with four distinct layers and a water resistant border. Comfortell combines an absorbent layer with a selectively permeable barrier that enables the wound to breathe while keeping out contaminants. Comfortell can be a primary dressing over a postoperative site, sutures or skin tears, and can also be applied as a secondary dressing with impregnated gauzes, wound fillers and enzymatic ointments.

Directions

1. Apply Gentell's Comfortell as a cover dressing for any primary or secondary treatment.
2. Flush the wound with Gentell Wound Cleanser and gently dry the skin surrounding the wound site.
3. As a primary dressing, apply directly to the wound surface.
4. As a secondary dressing, apply directly over primary treatment.
5. Repeat daily as ordered.