

## Your new Gentell Fastcare Tablet

Your account on [www.gentellfastcare.com](http://www.gentellfastcare.com) is

**LOGIN: XXXXXXXXXXXX**

**PASSWORD: XXXXXXXXXXXX**

Gentell has installed everything you need on your Gentell Fastcare tablet. *Installing additional software may interfere with its operation.*

You are using a **Motorola Xoom** tablet. Turn on the tablet by pressing the button on the back of the tablet behind the Motorola logo in the top left. Sound volume is also controlled by two buttons on the left side of the tablet.

When the screen turns on, **drag the picture of the lock** in any direction to unlock the screen. *When unused, the tablet turns off automatically after 5 minutes to conserve power.*

Tap on the icon of the **house** at the bottom left to go to the **Home screen** at any time.

Tap on the **Settings** icon, then **Wireless & networks**, then **Wi-Fi settings** to connect to your local Wi-Fi network. You may need a password to connect.

To access Fastcare (which requires Wi-Fi), tap on the **Firefox** icon on the home screen. If you are not already on the login screen, hold your finger on the top line until a list of web pages appears: choose the **Gentell Fastcare login** page (<http://www.gentellfastcare.com>).

*Because of HIPAA requirements, reps would not normally see information that is available to nurses.* The **demonstration account** (<http://demo.gentellfastcare.com>) enables Gentell reps to work like a nurse. Unlike the REAL Fastcare site, the demo account is not a secure web site, so you will see a warning when you first access the demo site.

Access this account with:

Email: **XXXXXXXXXXXX**

Password: **XXXXXXXXXXXX** (that's right: a hyphen and 2 t's)

If you'd like to play with the camera operations, please use the patient record called **Test Test**.

We suggest you watch the **Fastcare video demonstration** of the tablet on **gentell.com** for further instruction.

## **Teleconferencing**

You can use this tablet to teleconference, too. Tap on **Talk** on the home screen. This tablet's gtalk account is: **XXXXXXXXXXXX** with the password **XXXXXXXXXXXX**. If you have a web camera and a gtalk account on another computer, you can invite yourself or someone in your office to conference with you. When you are teleconferencing, you can tap the video camera icon to change to the camera that is pointing away from you; a nurse would use this camera when she is displaying a wound during a teleconference.

## **Showing Gentell documents on the tablet**

Tap on the **Adobe Reader** from the **Home screen**, and then tap on **PDF files** at the top right. You will see many documents that are helpful when speaking to customers about Gentell.

## **Shutting Down**

Tap the power button on the upper left for a quick shut off and turn on. To completely turn off the device (which requires additional time during the next start-up, but uses less battery power), hold down the power button for a moment until the screen goes black. If you do a total shutdown, you may need to start **sl4a** and the **send\_pics** script the next time you use Fastcare if you want take photos.