



Gentell Calcium Alginates are the best in class. Available in 4"x4", 6"x6" and 12" rope, Gentell Calcium Alginates were judged "**exceptional**" by an independent testing lab.

Qualifying facilities may trial Gentell Calcium Alginates at no charge. Please contact your Gentell Representative or call Gentell at 215-788-2700 to arrange a free trial.

The report by *SGS Life Science Services* is attached for your reference.



# CERTIFICATE OF ANALYSIS

Client: STEVE LEE  
GENTELL  
2701 BARTRAM ROAD  
BRISTOL PA 19007

Assignment No.: 09A1202847  
PO No.: VARFQ  
Page: 1 of 2

**Product:** Misc. Product

**Date Received:** 08-Mar-2012

**Qty Rec'd Per Sample No.:** 12

**Test/Method:** Analysis

SGS Sample No.	Identification	Method	Test Performed	Result	Status
09S12011768	4x4 Calcium Alginate Lot: 201112001 Sample Code: GEN-13500	BP Addendum 1995, page 1706	Absorbency	Exceptional	N/A

**Product:** Misc. Product

**Date Received:** 08-Mar-2012

**Qty Rec'd Per Sample No.:** 12

**Test/Method:** Analysis

SGS Sample No.	Identification	Method	Test Performed	Result	Status
09S12011769	12" Calcium Alginate Rope Lot: 201108001 GEN-13120	BP Addendum 1995, page 1706	Absorbency	Exceptional	N/A

This report is issued by SGS Life Science Services under its General Conditions for Testing Services (copy available upon request). SGS' responsibility under this report is limited to proven negligence and will in no case be more than the amount of the testing fees. Except by special arrangement, samples are not retained by SGS for more than 30 days. The results shown on this test report refer only to the sample(s) tested unless otherwise stated, under the conditions agreed upon. Anyone relying on this report should understand all of the details of the engagement. Neither the name, seals, marks nor insignia of SGS may be used in any advertising or promotional materials without the prior written approval of SGS. The test report cannot be reproduced, except in full, without prior written permission of SGS.



## CERTIFICATE OF ANALYSIS

Client: STEVE LEE  
GENTELL

Assignment No.: 09A1202847

Page: 2 of 2

### Evaluation / Comments

Acceptance criteria:

- Exceptional: Meets absorbency value of >11 g of solution per 100 cm<sup>2</sup> of sample.
- Good: Meets absorbency value of 10 to <11 g per 100 cm<sup>2</sup> of sample.
- Fair: Meets absorbency value of 9 to <10 g per 100 cm<sup>2</sup> of sample.
- Poor: Meets absorbency value of 8 to <9 g per 100 cm<sup>2</sup> of sample.

*This document was electronically approved by:*

*Theresa Larosiliere*

Date/Time: 12-Mar-2012, 03:40:36 PM