

**Advanced and  
Traditional  
Wound Care  
Product  
Reference  
Guide**



## Tables of Contents

---

<b>Super Absorbent Dressings</b>	Xtrasorb® Super Absorbent Dressings	<b>1</b>
	Gentell Super Absorbent Dressings	<b>2</b>
	Gentell Super Absorbent Adherent Dressings	<b>2</b>
<b>Collagen and Algicell Dressings</b>	Gentell Collagen	<b>3</b>
	Algicell® Ag Wound Dressings	<b>3</b>
<b>Calcium Alginate Dressing</b>	Algicell® Calcium Alginate Dressings	<b>4</b>
	Calcium Alginate Dressings	<b>4</b>
	Gentell Calcium Alginate Dressings	<b>5</b>
	Gentell Calcium Alginate Ag (Silver) Dressings	<b>5</b>
<b>CMC Fiber Dressings</b>	Gentell Fiber Dressings	<b>6</b>
<b>Wound Dressings</b>	Dermagran® B Hydrophilic Wound Dressings	<b>7</b>
	Dermagran® Hydrogel Wound Dressings	<b>7</b>
<b>Hydrogel Dressings</b>	Gentell Hydrogel Wound Dressings	<b>8</b>
	Aquasite Amorphous Hydrogel Dressings	<b>8</b>
	Aquasite Impregnated Gauzes	<b>8</b>
<b>Foam Dressing</b>	Sorbacell Foam Dressing (Without Film Backing)	<b>9</b>
	Gentell Waterproof Foam Dressing (Non-Bordered)	<b>9</b>
	Gentell Split Drain Foam Dressing (Non-Bordered)	<b>9</b>
	Gentell Bordered Foam	<b>10</b>
	Gentell Bordered Foam Ag (Silver)	<b>10</b>
	Gentell Bordered Silicone Foam	<b>10</b>
	Gentell Bordered Silicone Foam - Sacral	<b>11</b>
	Gentell Circular Split Drain Foam Dressing	<b>11</b>
<b>Hydrocolloid Dressings</b>	Dermatell Hydrocolloid Dressing (Non-Bordered)	<b>12</b>
	Dermatell Hydrocolloid Bordered Dressing	<b>12</b>
	Dermatell Hydrocolloid Dressing - Sacral	<b>12</b>
<b>Wound Cleansers</b>	Primaderm Wound Cleanser - Non-Sterile	<b>13</b>
	Primaderm Dermal Wound Cleanser - Sterile	<b>13</b>
	Dermagran® Wound Cleanser with Zinc	<b>13</b>
	Gentell Dermal Wound Cleanser	<b>14</b>
<b>Barrier Creams</b>	Shield & Protect™ Barrier Cream	<b>14</b>
	SuperMax Barrier Cream	<b>14</b>
<b>Sponges</b>	Ducare Woven Gauze Sponges	<b>15</b>
	Dupaque X-Ray Detectable Woven Gauze Sponge	<b>15</b>
	Dusoft Non-Woven Sponges	<b>16</b>
	Tracheotomy/Drain Sponges - Sterile	<b>16</b>

---

# Tables of Contents

---

<b>Wound Dressings</b>	Gazetex Bandage Rolls	<b>17</b>
	Gazetex 100% Cotton,6-Ply, Fluff Sponge	<b>17</b>
	Dupad Abdominal Pad	<b>18</b>
	Primapad Non-Adherent Absorbent Pad - Sterile	<b>18</b>
	Surecare Eye Pads - Sterile	<b>18</b>
<b>Packing</b>	Pak-Its Gauze Packing Strips - Plain and Iodoform	<b>19</b>
	Pak-Its Gauze Packing Strips - Plain and Iodoform with Strip Delivery System	<b>19</b>
<b>Island &amp; Barrier Border Gauze Dressings</b>	Compdress Island Dressing - Sterile	<b>20</b>
	Gentell Bordered Gauze	<b>20</b>
	Gentell Split Drain Bordered Gauze	<b>20</b>
	Gentell Rolled Gauze	<b>21</b>
	Dudress Film Top Barrier Dressing - Sterile	<b>21</b>
<b>Conforming Bandages</b>	Dutex 100% Cotton 2 Ply Conforming Bandage	<b>22</b>
	Duform Synthetic Conforming Bandage, Knitted	<b>22</b>
	Gentell Conforming Stretch Gauze	<b>23</b>
	Gentell Petrolatum Gauze	<b>24</b>
<b>Non-Adhering Impregnated Dressing</b>	Petrolatum Gauze Non-Adhering Dressing - Sterile	<b>24</b>
	Gentell Xeroform Dressing	<b>25</b>
	Xeroform Petrolatum Dressing - Sterile	<b>25</b>
	Oil Emulsion Dressing	<b>25</b>
	Shur-Conform Oil Emulsion Non-Adhering Dressing - Sterile	<b>26</b>
	MVP Transparent Film Dressing	<b>26</b>
	Comfortell Composite Dressing	<b>27</b>
	Surgilast <sup>®</sup> Tubular Elastic Dressing Retainer	<b>28</b>
<b>Composite Dressings</b>	Surgilast <sup>®</sup> Tubular Elastic Dressing Retainer	<b>29</b>
	Surgilast <sup>®</sup> Tubular Elastic Dressing Retainer	<b>30</b>
	Surgitube <sup>®</sup> Tubular Gauze (For use without applicator)	<b>31</b>
	Surgitube <sup>®</sup> Tubular Gauze (For use without applicator)	<b>32</b>
	SurgiGrip Tubular Elastic Support Bandage	<b>32</b>
	SurgiGrip Tubular Elastic Support Bandage	<b>33</b>
	Surgitube and SurgiGrip Tubular Gauze Applicators	<b>33</b>
<b>Elastic Net Dressing (Retainer/Tube Gauze)</b>		

---

## Tables of Contents

---

<b>Tapes/ Laparotomy Sponges</b>	Gentell FixTape	<b>34</b>
	Ultrafix Self-Adhesive Dressing Retention Tape	<b>34</b>
	Laparotomy Sponges, Prewashed and X-Ray Detectable	<b>34</b>
	Sorban Undercast Padding Bandage	<b>35</b>
<b>Compression Bandages &amp; Bandaging Systems</b>	Dufore Four-Layer Bandaging System	<b>35</b>
	Primer <sup>®</sup> Modified Unna Boot	<b>35</b>
	Unna-Pak: Primer <sup>®</sup> Modified Unna Boot and Duban Self Adherent Bandage	<b>36</b>
<b>Protective Products</b>	Glen-Sleeve II Arm and Leg Protectors	<b>37</b>
	Glen-Sleeve Arm Protectors	<b>37</b>
	Heel & Elbow Protectors	<b>37</b>
<b>Gauze Burn Dressings</b>	Dupress Burn Dressings	<b>38</b>
<b>Fixation Devices</b>	Cath-Strip <sup>®</sup> Reclosable Catheter Fastener	<b>39</b>
	UC Strip Catheter Tubing Fastener	<b>39</b>
	NG Strip <sup>®</sup> Nasal Tube Fastener	<b>39</b>
<b>Fixation Devices</b>	The Marsupial Pouch Wound Drainage Container Hol	<b>40</b>
<b>Wound Closure Strips</b>	Suture Strip <sup>®</sup> Plus Flexible Wound Closure Strips - Sterile	<b>41</b>
	Suture Strip <sup>®</sup> Wound Closure Strips - Non-Sterile	<b>41</b>
	Shur Strip Wound Closure Strips - Sterile	<b>42</b>
	Shur Strip Wound Closure Strips Non-Sterile Bulk	<b>42</b>

---

## Super Absorbent Dressings

### Xtrasorb<sup>®</sup> Super Absorbent Dressings

Xtrasorb<sup>®</sup> dressings come in a variety of formats and sizes. Each absorbs and retains more moisture than market leading competitive dressings. The exceptional fluid retention and moisture vapor transmission may contribute to increased time between dressing changes and reduce risk of maceration and exposure to the harmful components within wound fluid.



	SKU	Size	Packaging Unit/Case	HCPCS
HCS Non-Adhesive	86322	2.3" x 2.3"	10/box, 4 boxes/case	A6234
	86344	4.5" x 4.5"	10/box, 4 boxes/case	A6235
	86388	8" x 8"	5/box, 8 boxes/case	A6236
HCS Adhesive	86433	3" x 3"	10/box, 4 boxes/case	A6237
	86466	6" x 6"	10/box, 4 boxes/case	A6238
	86532	2.75" x 3.15"	10/box, 4 boxes/case	A6209
Gentle Tack Silicone Foam Borderless	86544	4.3" x 4.3"	10/box, 4 boxes/case	A6210
	86586	8.6" x 6"	5/box, 8 boxes/case	A6211
	86633	3.15" x 3.15"	10/box, 4 boxes/case	A6212
Gentle Tack Silicone Foam Bordered	86644	4" x 4"	10/box, 4 boxes/case	A6212
	86666	6" x 6"	10/box, 4 boxes/case	A6212
	86122	2" x 2"	10/box, 4 boxes/case	A6209
Foam Non-Adhesive	86144	4" x 4.75"	10/box, 4 boxes/case	A6210
	86188	8" x 8"	5/box, 8 boxes/case	A6211
	86233	3.2" x 3.2"	10/box, 4 boxes/case	A6212
Foam Adhesive	86244	4.5" x 4.5"	10/box, 4 boxes/case	A6212
	86266	6" x 6"	10/box, 4 boxes/case	A6213
	89533	3" x 3"	10/box, 10 boxes/case	A6251
Classic Non-Adhesive	89544	4" x 4"	10/box, 10 boxes/case	A6251
	89545	4" x 5"	10/box, 10 boxes/case	A6252
	89569	6" x 9"	10/box, 5 boxes/case	A6253

## Super Absorbent Dressings



### Gentell Super Absorbent Dressings

Gentell Super Absorbent Dressings offer excellent absorbent capacity for the treatment of moderate or heavy exuding wounds. The dressings consist of a breathable and waterproof non-woven backing, a super absorbent pad layer and a hydro-penetrating non-woven contact layer. Ideal for leg ulcers, pressure injuries, non-infected diabetic foot ulcers, dehisced surgical wounds and donor sites.

Super Absorbent  
Non-Adherent

SKU	Size	Packaging Unit/Case	HCPCS
GEN-19320	2" x 2"	10/box, 5 boxes/case	A6196/A6251
GEN-19340	4" x 4"	10/box, 5 boxes/case	A6196/A6251
GEN-19350	5" x 5"	10/box, 5 boxes/case	A6197/A6252
GEN-19380	8" x 8"	10/box	A6198/A6253



### Gentell Super Absorbent Adherent Dressings

Gentell Super Absorbent Adherent Dressings offer excellent absorbent capacity for the treatment of moderate or heavy exuding wounds. They provide all the benefits of Gentell's Super Absorbent Dressing with the addition of a comfortable silicone contact layer that provides gentle adherence and can speed up exudate absorption. The dressings consist of a breathable and waterproof non-woven backing, a super absorbent pad layer and a silicone contact layer. Ideal for leg ulcers, pressure injuries, non-infected diabetic foot ulcers, dehisced surgical wounds and donor sites.

Super Absorbent  
Adherent

SKU	Size	Packaging Unit/Case	HCPCS
GEN-19420	2" x 2"	10/box, 5 boxes/case	A6196
GEN-19450	5" x 5"	10/box, 5 boxes/case	A6197
GEN-19480	8" x 8"	10/box	A6198

## Collagen and Algicell Dressings

### Gentell Collagen

Gentell Collagen is a primary dressing for management of burns, scrapes, sores, ulcers, and blisters including wounds with minimal to heavy exudate, partial and full thickness, and second-degree burns.



	SKU	Size	Packaging Unit/Case	HCPCS
Collagen Dressings	GEN-18200	2" x 2"	10/box	A6021
	GEN-18400	4" x 4"	10/box	A6021
Collagen Particles	GEN-18700	7" x 7"	10/box	A6023
	GEN-18000	1g Pouch	5/box	A6010

### Algicell<sup>®</sup> Ag Wound Dressings

Algicell<sup>®</sup>Ag contains 1-4% ionic silver and provides a long-lasting, measured release for 7-day effectiveness. Algicell Ag absorbs 16%-32% more drainage than the market leader, reducing maceration risk. The dressings are effective against a broad range of bacteria, and the gelling capability ensures that the dressing doesn't stick or leave fibrous residue in the wound bed.



SKU	Size	Packaging Unit/Case	HCPCS
88512	.75" x 12"	5/box, 4 boxes/case	A6199
88522	2" x 2"	10/box, 10 boxes/case	A6196
88544	4.25" x 4.25"	10/box, 5 boxes/case	A6197
88545	4" x 5"	10/box, 5 boxes/case	A6197
88548	4" x 8"	5/box, 4 boxes/case	A6197
88582	8" x 12"	5/box, 4 boxes/case	A6198

## Calcium Alginate Dressings

### Algicell<sup>®</sup> Calcium Alginate Dressings

These soft, white sterile dressings absorb moderate to large amounts of wound drainage and may help control minor bleeding. Alginate ropes can be used to fill or pack wound. Algicell<sup>®</sup> dressings are designed to keep their form, allowing for easy and trauma-free removal, while maintaining a moist wound environment through its ability to gel upon contact with wound fluids. Does not contain natural rubber latex.



SKU	Size	Packaging Unit/Case	HCPCS
88022	2" x 2"	10/box, 4 boxes/case	A6196
88044	4" x 4"	10/box, 4 boxes/case	A6196
88048	4" x 8"	10/box, 4 boxes/case	A6197
88112	.75" x 12"	10/box, 4 boxes/case	A6199
88100	.75" x 36" Rope	1/box, 12 boxes/case	A6199

### Calcium Alginate Dressings

These soft, white sterile dressings absorb moderate to large amounts of wound drainage and may help control minor bleeding. Calcium alginate ropes can be used to fill or pack a wound. Calcium Alginate dressings are designed to keep their form, allowing for easy and trauma-free removal, while maintaining a moist wound environment through its ability to gel upon contact with wound fluids. Does not contain natural rubber latex.



SKU	Size	Packaging Unit/Case	HCPCS
88222	2" x 2"	10/box, 4 boxes/case	A6196
88222	4" x 4"	10/box, 4 boxes/case	A6196
88222	4" x 5"	10/box, 4 boxes/case	A6196
88222	4" x 8"	10/box, 4 boxes/case	A6196
88222	.75" x 12" Rope	10/box, 4 boxes/case	A6196



## Calcium Alginate Dressings



### Gentell Calcium Alginate Dressings

Gentell Calcium Alginate Dressings are a sterile, comfortable, advanced fiber-structured alginate with a highly absorbent capacity. Alginate dressings absorb, collect, and contain exudate while providing a moist healing environment. A reaction between the calcium in the dressing and the sodium in the wound exudate creates a gel-like substance that promotes moist wound healing. Use on leg ulcers, pressure ulcers, diabetic foot ulcers, and surgical wounds. May be used for minor conditions such as lacerations, abrasions, skin tears, and minor burns.

SKU	Size	Packaging Unit/Case	HCPCS
GEN-13200	2" x 2"	10/box, 5 boxes/case	A6196
GEN-13500	4" x 4"	10/box, 5 boxes/case	A6196
GEN-13480	4" x 8"	10/box, 10 boxes/case	A6197
GEN-13700	6" x 6"	10/box, 5 boxes/case	A6197
GEN-13010	10" x 10"	10/box, 5 boxes/case	A6198
GEN-13120	12" Rope	5/box, 10 boxes/case	A6199

### Gentell Calcium Alginate Ag (Silver) Dressings

Gentell Calcium Alginate Ag (Silver) Dressings are a sterile, antimicrobial, comfortable, fiber-structured alginate with high absorbency. These advanced alginate dressings fight bacteria and a broad spectrum of microorganisms while absorbing and containing exudate in a moist healing environment. A steady release of silver and a reaction between the calcium in the dressing and the sodium in the wound exudate creates a gel that fights infection and is easy to remove from the wound. May be used on wounds, such as leg ulcers, pressure ulcers, diabetic foot ulcers, and surgical wounds or minor conditions such as lacerations, abrasions, skin tears, and minor burns.



SKU	Size	Packaging Unit/Case	HCPCS
GEN-13220	2" x 2"	10/box, 10 boxes/case	A6196
GEN-13220	4" x 4"	10/box, 10 boxes/case	A6196
GEN-13220	4" x 8"	10/box, 10 boxes/case	A6197
GEN-13220	6" x 6"	10/box, 10 boxes/case	A6197
GEN-13220	8" x 8"	10/box, 10 boxes/case	A6198
GEN-13220	12" Rope	10/box, 10 boxes/case	A6199

## CMC Fiber Dressings



### Gentell Fiber Dressings

Gentell Fiber Dressing is a carboxymethylcellulose absorptive dressing for wounds with moderate to heavy exudate. A reaction between the fiber in the dressing and the wound fluid creates a gel-like substance that promotes moist wound healing and the formation of granulation tissue. The gelling fiber facilitates the removal of dead, damaged and infected tissue from the wound, and manages drainage to reduce the risk of periwound maceration.

SKU	Size	Packaging Unit/Case	HCPCS
GEN-13020	2" x 2"	10/box	A6196
GEN-13040	4" x 4"	10/box	A6196
GEN-13060	6" x 6"	10/box	A6197

## Wound Dressings



### Dermagran<sup>®</sup> B Hydrophilic Wound Dressings

These dressings provide either a primary cover or filler for chronic and acute wounds. Contains Derma Sciences' zinc-nutrient dressing formulation, balanced pH technology, vitamins A and B<sub>6</sub>, calcium, and magnesium. Indicated for skin ulcers (diabetic, venous stasis), stage I-IV pressure sores, surgical incisions and superficial injuries such as minor, partial thickness thermal burns, superficial lacerations, cuts or abrasions. Does not contain natural rubber latex.

SKU	Size	Packaging Unit/Case	HCPCS
SPD-03	3 oz Tube	12/case	A6261
SPD-20	2" x 2"	25/box, 8 boxes/case	A6231
SPD-21	4" x 4"	15/box, 4 boxes/case	A6231
SPD-24	4" x 8"	15/box, 4 boxes/case	A6232

### Dermagran<sup>®</sup> Hydrogel Wound Dressings



Hydrogel Wound Dressing, contains zinc and vitamin B<sub>6</sub>, may be used as a primary cover or filler for wound deficiencies. It absorbs wound exudate and creates a moist, mildly acidic environment to facilitate wound healing. Indicated for stage IV pressure ulcers, venous stasis ulcers, partial thickness thermal burns, surgical incisions, skin irritations, abrasions, cuts and skin conditions associated with peristomal care. Does not contain natural rubber latex.

SKU	Size	Packaging Unit/Case	HCPCS
DH-03	3 oz Tube	12/case	A6248

## Hydrogel Dressings



### Gentell Hydrogel Wound Dressings

Gentell Hydrogel is an Aloe Vera-based hydrating wound gel that protects the wound bed. Because it contains less water than other hydrogels, Gentell Hydrogel is more viscous and less “runny.” Aloe Vera is also a source of ace mannin, in addition to other mono- and polysaccharides, amino acids, glycoproteins, vitamins and enzymes. Use on wounds with little or no exudate, diabetic skin ulcers, venous stasis ulcers, first- and second-degree burns, post-surgical incisions, cuts, and abrasions.

SKU	Size	Packaging Unit/Case	HCPCS
GEN-11200	2" x 2"	10/box	A6231
GEN-11400	4" x 4"	10/box	A6231
GEN-11800	4" x 8"	10/box	A6232
GEN-11140	4 oz Tube	12/case	A6248
GEN-11080	8 oz Spray Bottle	6/case	N/A



### Aquasite Amorphous Hydrogel Dressings

Amorphous Hydrogels are useful for rehydrating sloughy or necrotic tissue on minimally draining wounds. They help to provide a moist wound environment, conducive to wound healing. Amorphous Hydrogel is used to fill uneven and difficult to reach spaces. The dressing is in a bellows tube, making precise application easier. Indicated for the management of acute and chronic wounds.

SKU	Size	Packaging Unit/Case	HCPCS
81110	1 fl oz Bellows	12/case	A6248




### Aquasite Impregnated Gauzes

Our impregnated dressings incorporate our Aquasite hydrogel into 100% Cotton USP Type VII gauze pads. The dressings are useful for rehydrating sloughy or necrotic tissue, helping to provide a moist wound environment on minimally draining wounds, and providing a sterile, moist packing for full thickness wounds. These dressings are ideal for the management of acute and chronic wounds including stage 4 pressure ulcers. Does not contain natural rubber latex.



SKU	Size	Packaging Unit/Case	HCPCS
53022	2" x 2"	10/box, 4 boxes/case	A6231
53044	4" x 4"	10/box, 4 boxes/case	A6231
53048	4" x 8"	10/box, 4 boxes/case	A6232

## Foam Dressing



### Sorbacell Foam Dressing (Without Film Backing)

Sorbacell<sup>®</sup> Foam dressing absorbs wound fluid and is flexible. A secondary dressing or medical tape is needed to hold the dressing in place. For use on moderately to heavily exuding wounds. Does not contain natural rubber latex.

SKU	Size	Packaging Unit/Case	HCPCS
84344	4" x 4"	10/box, 4 boxes/case	A6209



### Gentell Waterproof Foam Dressing (Non-Bordered)

Gentell Waterproof Non-Bordered Foam Dressing is a primary wound dressing for chronic and acute, moderate-to-heavy exuding, partial to full thickness wounds including superficial wounds and second-degree burns. Waterproof Non-Bordered Foam can also be used as a secondary wound dressing.

SKU	Size	Packaging Unit/Case	HCPCS
GEN-14522	2" x 2"	10/box, 5 boxes/case	A6209
GEN-14544	4" x 4"	10/box, 5 boxes/case	A6209
GEN-14555	5" x 5"	10/box, 5 boxes/case	A6210
GEN-14588	8" x 8"	10/box	A6211

### Gentell Split Drain Foam Dressing (Non-Bordered)

Gentell Split Drain Foam Dressing is a pre-cut, non-bordered, rectangular dressing with a U-shaped fenestration ideal for tracheostomy sites, catheters, feeding tubes or larger wound sites requiring moderate to heavy absorption. This highly absorptive foam dressing reduces strike-through and prevents skin maceration by keeping the wound site dry. May be cut and applied to pressure ulcers, chronic and acute wounds, and full thickness wounds.

SKU	Size	Packaging Unit/Case	HCPCS
GEN-14540	4" x 5"	10/box, 5 boxes/case	A6210



## Foam Dressing

### Gentell Bordered Foam

Gentell Bordered Foam Dressings are waterproof and highly absorbent foam dressings with an island design. Their dense, sturdy structure protects the wound from external threats, and the porous texture is highly absorbent. A nonabrasive border holds the foam in place and creates a water-resistant barrier to protect the wound. May be used as a primary or secondary dressing.



SKU	Size	Packaging Unit/Case	HCPCS
GEN-14400	4" x 4"	10/box, 5 boxes/case	A6212
GEN-14600	6" x 6"	10/box, 5 boxes/case	A6212

### Gentell Bordered Foam Ag (Silver)

Gentell Bordered Foam with Silver (Ag) is a silver ion-containing foam dressing with excellent anti-bacterial effect. It is ideal for the treatment of moderate or heavy exudate chronic and acute wounds caused by bacterial infection, including leg ulcers, pressure injuries, diabetic foot ulcers, second degree burns, donor sites, post-operative wounds and skin abrasions.



SKU	Size	Packaging Unit/Case	HCPCS
GEN-14844	4" x 4"	10/box, 5 boxes/case	A6212
GEN-14866	6" x 6"	10/box, 5 boxes/case	A6212

### Gentell Bordered Silicone Foam

Gentell Silicone Foam Bordered Dressing is a primary or secondary wound dressing for chronic and acute, partial, and full thickness wounds including superficial wounds and second-degree burns. May be used as a primary or secondary dressing.



SKU	Size	Packaging Unit/Case	HCPCS
GEN-14744	4" x 4"	10/box, 5 boxes/case	A6212
GEN-14766	6" x 6"	10/box, 5 boxes/case	A6212
GEN-14777	7" x 7"	10/box, 5 boxes/case	A6213
GEN-14799	9" x 9"	10/box, 5 boxes/case	A6214

## Foam Dressing



### Gentell Bordered Silicone Foam - Sacral

Gentell Silicone Foam Bordered Sacral Dressing is a primary or secondary wound dressing for chronic and acute, partial, and full thickness wounds including superficial wounds and second-degree burns. May be used as a primary or secondary dressing.

SKU	Size	Packaging Unit/Case	HCPCS
GEN-14700	9" x 9"	10/box	A6213



### Gentell Circular Split Drain Foam Dressing

Gentell Circular Split Drain Foam Dressing is a pre-cut circular, bordered, absorptive dressing with a U-shaped fenestration ideal for heels, joints, and other wound sites with irregular or protruding surfaces, plus ostomy sites, catheters, and feeding tubes. This super-absorptive foam dressing reduces strike-through and prevents skin maceration by keeping the wound site dry. The hypo-allergenic, water-resistant tape allows the skin to breathe and makes application easy, which reduces nursing time.

SKU	Size	Packaging Unit/Case	HCPCS
GEN-14410	4"	10/box, 5 boxes/case	A6212

# Hydrocolloid Dressings

## Dermatell Hydrocolloid Dressing (Non-Bordered)

Gentell Dermatell Hydrocolloid Dressings consist of a soft, pliable hydrocolloid wafer that enhances patient comfort and protection. Dermatell is most effective when kept in place for a minimum of three days. Dermatell is available with a water-resistant adhesive border that is flexible and conforms easily to the body, and in a non-bordered hydrocolloid that naturally adheres to wounds with exudate.



SKU	Size	Packaging Unit/Case	HCPCS
GEN-10100	4" x 4"	10/box, 5 boxes/case	A6234

## Dermatell Hydrocolloid Bordered Dressing

SKU	Size	Packaging Unit/Case	HCPCS
GEN-10200	4" x 4"	10/box, 5 boxes/case	A6237
GEN-10400	6" x 6"	10/box, 5 boxes/case	A6238



## Dermatell Hydrocolloid Dressing - Sacral

SKU	Size	Packaging Unit/Case	HCPCS
GEN-11100	6" x 7"	10/box, 5 boxes/case	A6238





## Wound Cleaners



### Primaderm Wound Cleanser - Non-Sterile

Non-sterile, preserved dermal cleanser for multiple use. Non-cytotoxic formula. Use when extra cleansing power is needed. Contains Poloxamer 188, an effective surfactant known to be gentle on tissue. Indicated for pressure ulcers, stasis ulcers, diabetic foot ulcers, post-surgical wounds, first and second degree burns, cuts, abrasions, and minor irritations of the skin.

SKU	Size	Packaging Unit/Case	HCPSC
69201	16 fl oz Spray Bottle	12/case	A6260
69202	4.15 oz Squirt Top Bottle	12/case	A6260



### Primaderm Dermal Wound Cleanser - Sterile

Sterile, non-cytotoxic dermal cleanser for single-use, making it the cleanser of choice for traumatic acute wounds. Use when extra cleansing power is needed. Contains Poloxamer 188, an effective surfactant known to be gentle on tissue. May be used for pressure ulcers, stasis ulcers, diabetic foot ulcers, post-surgical wounds, first and second degree burns, cuts, abrasions, and minor irritations of the skin.

SKU	Size	Packaging Unit/Case	HCPSC
69101	4.15 oz Squirt Top Bottle	12/case	A6260



### Dermagran<sup>®</sup> Wound Cleanser with Zinc

Dermagran<sup>®</sup> Wound Cleanser is a gentle skin and wound cleanser that contains both Zinc and Vitamin B6. This wound cleanser is designed to not just cleanse, but to also prepare the wound bed as an environment conducive to wound healing. The multi-use spray bottle is designed so that the protective cap does not need to be fully removed in order to operate, increasing ease of use. This 360° delivery system meets AHCPH guidelines for irrigation pressure.

SKU	Size	Packaging Unit/Case	HCPSC
WC-04	4 oz Bottle	12/case	A6260

## Wound Cleansers



### Gentell Dermal Wound Cleanser

Gentell Dermal Wound Cleanser is a no rinse, non-irritating cleanser for efficient, thorough, and gentle cleansing of the skin. Dermal Wound Cleanser can be used on any wound, regardless of stage.

SKU	Size	Packaging Unit/Case
GEN-10080	8 oz Spray Bottle	6/case
GEN-10160	16 oz Spray Bottle	6/case



### Shield & Protect™ Barrier Cream

Gentell Shield & Protect Moisture Barrier Cream forms a water-shedding shield that helps soothe and protect dry or irritated skin. Therapeutic against irritation from urine, feces, and other bodily secretions. Shield & Protect leaves skin soft and smooth rather than irritated and greasy. Ideal for sensitive skin.

SKU	Size	Packaging Unit/Case
GEN-23140	4 oz Tube	12/case



### SuperMax Barrier Cream

Gentell SuperMax Moisture Barrier Cream protects healthy and intact skin from breaking down during compromised skin conditions. SuperMax locks out moisture with 16% Zinc Oxide and Balsam of Peru by forming a water-resistant barrier to soothe and protect irritated skin. Contains a gentle, natural vanilla scent.

SKU	Size	Packaging Unit/Case
GEN-23340	4 oz Tube	12/case

## Sponges

### Ducare Woven Gauze Sponges



These woven sponges are made from 100% cotton and may be used for general use, covering and packing wounds. These dressings come in a variety of sizes and are available as sterile or non-sterile. Does not contain natural rubber latex.

	SKU	Size	Packaging Unit/Case
Sterile	82208	2" x 2", 8 ply	2 ea/pouch, 50 pouches/box, 30 boxes/case, 3000 ea/case
	82312	3" x 3". 12 ply	2 ea/pouch, 25 pouches/box, 48 boxes/case, 2400 ea/case
	82408	4" x 4", 8 ply	2 ea/pouch, 25 pouches/box, 24 boxes/case, 1200/case
	82410	4" x 4", 12 ply	10 ea/pouch, 32 pouches/tray. 4 trays/case, 1280 ea/case
	82412	4" x 4", 12 ply	2 ea/pouch, 25 pouches/box. 24 boxes/case, 1200 ea/case
	82812	8" x 4", 12 ply	2 ea/pouch, 25 pouches/box. 20 boxes/case, 1000 ea/case
Non-sterile	90208	2" x 2", 8 ply	200/sleeve, 25 sleeves/case. 5000/case
	90212	2" x 2", 12 ply	200/sleeve, 40 sleeves/case. 8000/case
	90312	3" x 3", 12 ply	200/sleeve, 20 sleeves/case. 4000/case
	90408	4" x 4", 8 ply	200/sleeve, 20 sleeves/case. 4000/case
	90412	4" x 4", 12 ply	200/sleeve, 10 sleeves/case, 2000/case
	90812	8" X 4", 12 ply (H Fold)	200/sleeve, 10 sleeves/case 2000/case

### Dupaque X-Ray Detectable Woven Gauze Sponge



These woven sponges are made from 100% cotton and are low-linting, fluffy, and highly absorbent. They feature x-ray detectable barium sulfate strands for easy detection. Does not contain natural rubber latex.

	SKU	Size	Packaging Unit/Case
Sterile	9464S	4" x 4", 16 ply	10 ea/tray, 128 trays/case, 1280/case
	9467S	8" x 4". 12 ply	10 ea/tray, 80 trays/case. 800/case

## Sponges



### Dusoft Non-Woven Sponges

These non-woven sponges are ideal for wound cleansing and also can be used as a secondary dressing for added absorbency. They produce less lint and are less adherent than woven gauze sponges, so they cause less disruption to the wound bed. The rayon/polyester blend provides fast wicking properties while being highly absorbent, and minimizing skin maceration. Does not contain natural rubber latex.

	SKU	Size	Packaging Unit/Case
Sterile	84122	2" x 2", 4 ply	2 ea/pouch, 25 pouches/box, 60 boxes/case, 3000ea/case
	84133	3" x 3", 4 ply	2 ea/pouch, 25 pouches/box, 48 boxes/case, 2400 ea/case
	84144	4" x 4", 4 ply	2 ea/pouch, 25 pouches/box, 24 boxes/case, 1200ea/case
	84644	4" x 4", 6 ply	2 ea/pouch, 25 pouches/box, 12 boxes/case, 600 ea/case
Non-sterile	94122	2" x 2", 4 ply	200 ea/sleeve, 40 sleeves/case, 8000 ea/case
	94133	3" x 3", 4 ply	200 ea/sleeve, 20 sleeves/case, 4000 ea/case
	94144	4"x 4", 4 ply	200 ea/sleeve, 10sleeves/case, 2000 ea/case



### Tracheotomy/Drain Sponges - Sterile

Added fenestration makes these sponges ideal for dressing around a tracheostomy tube or any type of tube on the body. Available in a soft 6 ply non-woven or in standard 12 ply gauze. These sponges are absorbent, low-linting, fast-wicking, and non-adherent to wound beds. The woven sponges are ideal for when natural breathability is desired. Does not contain natural rubber latex.

SKU	Size	Packaging Unit/Case
84916	Non-Woven 4" x 4", 6-Ply	2's, 25 pkgs/box, 12 boxes/case, 600/case
84918	Non-Woven 2" x 2", 6-Ply	2's, 35 pkgs/box, 20 boxes/case, 1400/case
82910	Woven Gauze 4" x 4", 12-Ply	1's 25 pkgs/box, 24 boxes/case, 600/case

## Wound Dressings

### Gazetex Bandage Rolls

Gazetex bandage rolls are washed, low-linting woven gauze rolls. These bandages are used for wrapping or packing large and deep wounds. Does not contain natural rubber latex. being highly absorbent, and minimizing skin maceration. Does not contain natural rubber latex.



	SKU	Size	Packaging Unit/Case
Sterile	9322	4-1/2"x 4.1yds, 6 ply	1/pkg, 100 pkgs/case, 100/case
	9321	2-1/2"x 3 yds, 6 ply	1/pkg, 100 pkgs/case, 100/case
Non-sterile	9320	4-1/2"x 4.1yds, 6 ply	1/pkg, 100 pkgs/case, 100/case
	9320-48	4-1/2" x 4.1yds, 6 ply	12/sleeve, 4 sleeves/case, 48/case

### Gazetex 100% Cotton,6-Ply, Fluff Sponge

These washed, fluffed gauze sponges are low-linting and soft. They are excellent for absorption, packing wounds, and managing exudate in and around the wound. The dressings are ideal for pre-op prepping and cleansing. Does not contain natural rubber latex.



	SKU	Size	Packaging Unit/Case
Sterile	9340	6" x 6l-3/4", diagonally	2's, 240 pkgs/case, 480/case
	9345	7-3/4"x 8-3/4", diagonally	2's, 80 pkgs/case, 160/case
Non-sterile	9333	6" x 6-3/4", diagonally	100/sleeve, 6 sleeves/case, 600/case
	9334	7-3/4" x 8-3/4", diagonally	100/sleeve, 4 sleeves/case, 400/case

## Wound Dressings



### Dupad Abdominal Pad

Our abdominal pads are filled with cellulose fibers, making them super-absorbent. The fibers quickly absorb and disperse fluids laterally, to prevent pooling and maceration. The backing is constructed of a hydrophobic material to prevent strike-through. All edges are sealed to prevent linting. Does not contain natural rubber latex

	SKU	Size	Packaging Unit/Case
Sterile	87059	5" x 9"	25/box, 8 boxes/case, 200/case
	87088	8" x 7-1/2"	25/box, 8 boxes/case, 200/case
	87810	8"x10"	25/box, 8 boxes/case, 200/case
Non-sterile	97159	5" x 9"	576/case (Bulk packaged)
	97188	8" x 7-1/2"	576/case (Bulk packaged)
	97110	8" x 10"	432/case (Bulk packaged)



### Primapad Non-Adherent Absorbent Pad - Sterile

This non-adherent gauze pad is often chosen as the primary dressing when extra wound bed protection is needed. Air circulates through the soft, perforated film, while fluid is absorbed by the soft cotton pad. These dressings can be easily cut to fit any shape. Does not contain natural rubber latex.

SKU	Size	Packaging Unit/Case
83323	2" 3"	100/box, 24 boxes/case, 2400/case
83334	3" x 4"	50/box, 18 boxes/case, 900/case
83338	3" x 8"	50/box, 12 boxes/case, 600/case



### Surecare Eye Pads - Sterile

Constructed of absorbent cotton filler covered with fine mesh gauze on both sides. This pad is ideal for applying medication or cover and protect the eye. Does not contain natural rubber latex.

SKU	Size	Packaging Unit/Case
82911	2 5/8" x 1 5/8"	1's, 25/tray, 24 trays/case, 600/case

## Packing



### Pak-Its Gauze Packing Strips - Plain and Iodoform

Plain and Iodoform Gauze Strips are sterile and may be used dry to absorb drainage in exuding wounds or soaked with a medication of choice. Pre-cut sizes suit different wound types. A secondary dressing is needed to hold the strips in place and to protect the wound from outside contaminants. Does not contain natural rubber latex.

Plain Packing Strips

Iodoform Packing Strips

SKU	Size	Packaging Unit/Case
59121	1/4"X 5 yd	12/case
59221	1/2"x 5 yd	12/case
59321	1" x 5 yd	12/case
59421	2" x 5 yd	12/case
59146	1/4"x 5 yd	12/case
59246	1/2"x 5 yd	12/case
59346	1" x 5yd	12/case



### Pak-Its Gauze Packing Strips - Plain and Iodoform with Strip Delivery System

These packing strips are identical to our other packing strips except the bottles each have a special strip delivery system that ensures packing material is kept clean and tangle free, with the end always easy to locate at the top of the bottle. Does not contain natural rubber latex.

Plain Packing Strips

Iodoform Packing Strips

SKU	Size	Packaging Unit/Case
59120	1/4"X 5 yd	12/case
59220	1/2"x 5 yd	12/case
59320	1" x 5 yd	12/case
59420	2" x 5 yd	12/case
59145	1/4"x 5 yd	12/case
59245	1/2"x 5 yd	12/case
59345	1" x 5 yd	12/case
59445	2" x 5 yd	12/case

## Island & Barrier Border Gauze Dressings



### Compress Island Dressing - Sterile

These border gauze dressings can be used as either primary or secondary dressing. The product is comprised of a non-adherent, absorptive central pad, with a soft conformable layer of gauze as a backing. There is a skin-friendly adhesive border around the perimeter. Does not contain natural rubber latex.

SKU	Size	Packaging Unit/Case
89022	2" x 2"	25/box, 4 boxes/case 100/case
89044	4" x 4"	25/box, 4 boxes/case 100/case
89066	6" x 6"	25/box, 4 boxes/case 100/case
89046	4" x 6"	25/box, 4 boxes/case 100/case
89048	4" x 8"	25/box, 4 boxes/case 100/case
89041	4" x 10"	25/box, 4 boxes/case 100/case
89068	6" x 8"	25/box, 4 boxes/case 100/case
89042	4" x 14"	25/box, 4 boxes/case 100/case



### Gentell Bordered Gauze

Gentell Bordered Gauze is a multi-purpose wound dressing that saves nursing time and is easy to apply. A non-adherent pad is paired with a conforming beveled spunlace tape to create a secure, comfortable dressing that stays in place. May be used as a primary or secondary dressing.

SKU	Size	Packaging Unit/Case	HCPSC
GEN-15410	4" x 4"	25/box, 4 boxes/case	A6219
GEN-15610	6" x 6"	25/box, 4 boxes/case	A6220
GEN-15810	8" x 8"	25/box, 4 boxes/case	A6221

### Gentell Split Drain Bordered Gauze

Gentell Circular Split Drain Gauze Dressing is pre-cut circular, bordered dressing with a U-shaped fenestration ideal for heels, joints, and other wound sites with irregular or protruding surfaces, plus ostomy sites, catheters and feeding tubes.



SKU	Size	Packaging Unit/Case	HCPSC
GEN-15440	4"	10/box, 5 boxes/case	A6219



## Island & Barrier Border Gauze Dressings

### Gentell Rolled Gauze

Gentell Rolled Gauze Bandage is a sterile, conformable, and stretchable gauze bandage roll. The weave construction allows for a flexible and easy to use dressing designed to suit a variety of wound dimensions and locations.

SKU	Size	Packaging Unit/Case	HCPCS
GEN-15450	4.5" x 4.1 yards	10/box	A6446



### Dudress Film Top Barrier Dressing - Sterile

These all-in-one dressings are identical to our Compdress Island Gauze Dressings with the addition of polyurethane film barrier. The film barrier keeps out external contaminants as well as being an effective water barrier to permit showering and bathing. Does not contain natural rubber latex.



SKU	Size	Packaging Unit/Case
89144	4" x 4"	25/box, 4 boxes/case 100/case
89166	6" x 6"	25/box, 4 boxes/case 100/case
89168	6" x 8"	25/box, 4 boxes/case 100/case

## Conforming Bandages



### Dutex 100% Cotton 2 Ply Conforming Bandage

This 100% cotton 2-ply woven bandage is ideal as a secondary dressing for wrapping legs and arms and to hold dressings in place. The dressing provides light stretch which enables it to conform to difficult body contours. While also providing light compression. Does not contain natural rubber latex.

	SKU	Size	Packaging Unit/Case
Sterile	77781	2" x 4.1 yd	12/sleeve, 8 sleeves/case, 96/case
	77782	3" x 4.1 yd	12/sleeve, 8 sleeves/case, 96/case
	77783	4" x 4.1 yd	12/sleeve, 8 sleeves/case, 96/case
	77784	6" x 4.1 yd	6/sleeve, 8 sleeves/case, 48/case
Non-sterile	76781	2" x 4.1 yd	12/sleeve, 8 sleeves/case, 96/case
	76782	3" x 4.1 yd	12/sleeve, 8 sleeves/case, 96/case
	76783	4" x 4.1 yd	12/sleeve, 8 sleeves/case, 96/case
	76784	6" x 4.1 yd	6/sleeve, 8 sleeves/case, 48/case



### Duform Synthetic Conforming Bandage, Knitted


This knitted stretch gauze roll is ideal as a secondary dressing for wrapping legs and arms and to hold dressings in place. The dressings gently secure to body contours without restricting movement. Does not contain natural rubber latex.

	SKU	Size	Packaging Unit/Case
Sterile	75202	2" x 4.1 yd stretched	12/sleeve, 8 sleeves/case, 96/case
	75203	3" x 4.1 yd stretched	12/sleeve, 8 sleeves/case, 96/case
	75204	4" x 4.1 yd stretched	12/sleeve, 8 sleeves/case, 96/case
	75206	6" x 4.1 yd stretched	6/sleeve, 8 sleeves/case, 48/case
Non-sterile	75102	2" x 4.1 yd stretched	12/sleeve, 8 sleeves/case, 96/case
	75103	3" x 4.1 yd stretched	12/sleeve, 8 sleeves/case, 96/case
	75104	4" x 4.1 yd stretched	12/sleeve, 8 sleeves/case, 96/case
	75106	6" x 4.1 yd stretched	6/sleeve, 8 sleeves/case, 48/case

## Conforming Bandages

### Gentell Conforming Stretch Gauze

Gentell Conforming Stretch Gauze is a sterile, stretchable, and easy to manipulate gauze bandage roll. The weave construction allows for a flexible and easy to use dressing designed to suit a variety of wound dimensions and locations.



SKU	Size	Packaging Unit/Case	HCPCS
GEN-15310	3" x 75"	30/box	A6446

## Non-Adhering Impregnated Dressing



### Gentell Petrolatum Gauze

Gentell Petrolatum Dressing is an absorbent fine mesh gauze impregnated with whole petrolatum. Latex free and non-adherent.

SKU	Size	Packaging Unit/Case	HCPCS
GEN-19100	3" x 9"	50/box	A6223

### Petrolatum Gauze Non-Adhering Dressing - Sterile

Petrolatum-impregnated dressings are non-adherent and provide a moist environment that is conducive to wound healing for full-thickness wounds. They are made of fine, soft, conformable gauze-impregnated with white petrolatum. Indicated for the management of stage II-IV pressure ulcers, tunneling wounds, and non-infected wounds. Does not contain natural rubber latex.



SKU	Size	Packaging Unit/Case
DKC20049	1/2" x 72"	12/box, 6 boxes/case 72/case
DKC20053	1" x 36"	12/box, 6 boxes/case 72/case
DKC20056	3" x 9"	12/box, 6 boxes/case 72/case
DKC20059	3" x 18"	12/box, 6 boxes/case 72/case
DKC20062	3"x 36"	12/box, 6 boxes/case 72/case
DKC20065	6" x 36"	12/box, 6 boxes/case 72/case
DKC77030	1" x 8"	50/box, 4 boxes/case 200/case

# Non-Adhering Impregnated Dressing



## Gentell Xeroform Dressing

Gentell Xeroform Dressing is an absorbent, fine mesh gauze impregnated with 3% Bismuth Tribromophenate in a petrolatum blend. Latex free and non-adherent.

SKU	Size	Packaging Unit/Case	HCPCS
GEN-19210	1" x 8"	25/box	A6222
GEN-19220	2" x 2"	25/box	A6222
GEN-19240	4" x 4"	50/box	A6222
GEN-19200	5" x 9"	50/box	A6223

## Xeroform Petrolatum Dressing - Sterile

Xeroform impregnated dressings are similar to our Petrolatum Gauze dressings, with the addition of 3% bismuth tribromophenate. These dressings are for when mild medication or deodorizing are required, such as post-operative applications, 1st and 2nd degree burns, and skin grafts. Does not contain natural rubber latex.



SKU	Size	Packaging Unit/Case
DKC77034	1" x 8"	50/box, 4 boxes/case, 200/case
DKC20068	5" x 9"	50/box, 4 boxes/case, 200/case
DKC71164	2" x 2"	25/box, 6 boxes/case, 150/case
DKC71165	4" x 4"	25/box, 6 boxes/case, 150/case


## Oil Emulsion Dressing

Gentell Oil Emulsion Dressing with USP mineral oil is a conformable dressing ideal for lightly exuding wounds including minor burns, lacerations, and abrasions. Latex free and non-adherent.



SKU	Size	Packaging Unit/Case	HCPCS
GEN-19000	3" x 3"	25/box	A6222

## Non-Adhering Impregnated Dressing



### Shur-Conform Oil Emulsion Non-Adhering Dressing - Sterile

Oil Emulsion impregnated dressings are non-adherent and provide a moist environment that is conducive to wound healing. These dressings are made of a knitted cellulose acetate fabric coated with a formulated petrolatum emulsion. The knitted fabric allows wound fluid to move through the dressing and into a secondary absorbent dressing. Indicated for the management of skin grafts, surgical sites, abrasions, lacerations, dermal ulcers, and cosmetic surgery sites. Does not contain natural rubber latex.

SKU	Size	Packaging Unit/Case
DKC77041	3" x 3"	50/box, 12 boxes/case, 600/case
DKC77042	3" x 8"	24/box, 6 boxes/case, 144/case
DKC77043	3" x 16"	36/box, 6 boxes/case, 216/case
DKC77044	3" x 8", 3's	36/box, 6 boxes/case, 216/case

### MVP Transparent Film Dressing

Gentell MVP Transparent Film Dressing is a moisture vapor permeable transparent membrane coated with a layer of acrylic, hypo-allergenic adhesive that can be used to cover low exuding wounds. Gentell MVP is a non-absorptive sterile dressing that is permeable to moisture vapor and oxygen, but impermeable to bacteria. It can be used as a secondary dressing for Stage 2 pressure ulcers, abrasions, skin tears, blisters, skin graft donor sites, superficial partial-thickness burns, autolytic debridement, skin protection against moisture as well as friction and clean, closed surgical incisions. MVP Transparent dressings may also be used to cover and secure I.V. devices.

SKU	Size	Packaging Unit/Case	HCPCS
GEN-16450	4" x 4.75"	50/box	A6258


REORDER  
GEN-16450

MVP DRESSING

Transparent Film 4" x 4.75" (10cm x 12cm)

*Non-absorptive Moisture Vapor Permeable Dressing with hypoallergenic adhesive to cover and monitor low exuding wounds*

**Directions**

1. Flush with wound with Gentell Wound Cleanser and gently dry the skin surrounding the wound.
2. As a primary dressing, apply directly to the wound surface.
3. As a secondary dressing, apply directly over primary treatment.
4. Change daily or as ordered by a physician.



2001 Berwyn Road, Berwyn, PA 19007  
717-363-3111 FAX: 717-363-3170  
www.gentell.com

16298

## Composite Dressings



### Comfortell Composite Dressing

Gentell Comfortell<sup>®</sup> Dressing is a composite wound dressing with four distinct layers and a water-resistant border. Comfortell combines an absorbent layer with a selectively permeable barrier that enables the wound to breathe while keeping out contaminants. Comfortell can be a primary dressing over a postoperative site, sutures or skin tears, and can also be applied as a secondary dressing with impregnated gauzes, wound fillers and enzymatic ointments. May be used as a primary or secondary dressing.

SKU	Size	Packaging Unit/Case	HCPCS
GEN-12400	4" x 4"	10/box, 5 boxes/case	A6203
GEN-12600	6" x 6"	10/box, 5 boxes/case	A6203

## Elastic Net Dressing (Retainer/Tube Gauze)



### Surgilast® Tubular Elastic Dressing Retainer

A tubular elastic stretch net is designed to serve as a secondary dressing. It applies gentle pressure to keep bandages or other medical devices securely in place without adhesive tape. Ideal for securing burn or post-op dressings, hot or cold packs, IV tubes, incontinence pads and decubitus ulcer dressings. Easy to customize. Allows maximum air-flow.

Product	Size	Working Stretch	Applications	Pkg Unit/ Case
GL501	1	6 7/8"	Fingers, toes, wrist	50 yards, 1 box
GL702	2	8"	Small: hand, arm, leg, foot	25 yards, 1 box
GL502	2	8"	Small: hand, arm, leg, foot	50 yards, 1 box
GL720	3	10 1/8"	Medium: hand, arm, leg, foot	10 yards, 1 box
GL703	3	10 1/8"	Medium: hand, arm, leg, foot	25 yards, 1 box
GL503	3	10 1/8"	Medium: hand, arm, leg, foot	50 yards, 1 box
GL704	4	12 1/4"	Large: hand, arm, leg, foot	25 yards, 1 box
GL504	4	12 1/4"	Large: hand, arm, leg, foot	50 yards, 1 box
GL722	5	15"	Small: head, shoulder, thigh	10 yards, 1 box
GL705	5	15"	Small: head, shoulder, thigh	25 yards, 1 box
GL505	5	15"	Small: head, shoulder, thigh	50 yards, 1/box
GL706	5-1/2	19 3/4"	Medium: head, shoulder, thigh	25 yards, 1/box
GL506	5-1/2	19 3/4"	Medium: head, shoulder, thigh	50 yards, 1/box
GL724	6	25 1/2"	Large: head, shoulder, thigh	10 yards, 1/box
GL707	6	25 1/2"	Large: head, shoulder, thigh	25 yards, 1/box
GL507	6	25 1/2"	Large: head, shoulder, thigh	50 yards, 1/box
GL708	7	29"	Small: chest, back, perineum, axilla	25 yards, 1/box
GL508	7	29"	Small: chest, back, perineum, axilla	50 yards, 1/box
GL709	8	32 1/2"	Medium: chest, back, perineum, axilla	25 yards, 1/box
GL509	8	32 1/2"	Medium: chest, back, perineum, axilla	50 yards, 1/box
GL710	9	36"	Large: chest, back, perineum, axilla	25 yards, 1/box
GL510	9	36"	Large: chest, back, perineum, axilla	50 yards, 1/box
GL726	10	38"	X-Large: chest, back, perineum, axilla	10 yards, 1/box
GL711	10	38"	X-Large: chest, back, perineum, axilla	25 yards, 1/box
GL511	10	38"	X-Large: chest, back, perineum, axilla	50 yards, 1/box
GL712	11	40"	Special Sizing	25 yards, 1/box

Continued on next page



## Elastic Net Dressing (Retainer/Tube Gauze)

### Surgilast<sup>®</sup> Tubular Elastic Dressing Retainer (continued)

Product	Size	Working Stretch	Applications	Pkg Unit/Case
GL512	11	40"	Special Sizing	50 yards, 1/box
GL713	12	42"	Special Sizing	25 yards, 1/box
GL714	13	44"	Special Sizing	25 yards, 1/box
GL715	14	46"	Special Sizing	25 yards, 1/box
GL2522	22	72"	Special Sizing	25 yards, 1/box

### Surgilast<sup>®</sup> Tubular Elastic Dressing Retainer

Provides the same features and benefits as Surgilast while eliminating the latex sensitivity. Does not contain natural rubber latex.



Product	Size	Working Stretch	Applications	Pkg Unit/Case
GLLF2501	1	5 3/8"	Fingers, toes, wrist	25 yards, 1/box
GLLF2502	2	7"	Small: hand, arm, leg, foot	25 yards, 1/box
GLLF2503	3	9"	Medium: hand, arm, leg, foot	25 yards, 1/box
GLLF2504	4	11"	Large: hand, arm, leg, foot	25 yards, 1/box
GLLF2505	5	13 3/4"	Small: head, shoulder, thigh	25 yards, 1/box
GLLF2505.5	5 1/2	18"	Medium: head, shoulder, thigh	25 yards, 1/box
GLLF2506	6	23"	Large: head, shoulder, thigh	25 yards, 1/box
GLLF2507	7	28"	Small: chest, back, perineum, axilla	25 yards, 1/box
GLLF2508	8	31 1/2"	Medium: chest, back, perineum, axilla	25 yards, 1/box
GLLF2509	9	33"	Large: chest, back, perineum, axilla	25 yards, 1/box
GLLF2510	10	37 1/2"	X-Large: chest, back, perineum, axilla	25 yards, 1/box
GLLF2511	11	40"	Special Sizing	25 yards, 1/box
GLLF2512	12	41 1/2"	Special Sizing	25 yards, 1/box
GLLF2514	14	45"	Special Sizing	25 yards, 1/box

## Elastic Net Dressing (Retainer/Tube Gauze)



### Surgilast<sup>®</sup> Tubular Elastic Dressing Retainer

A pre-cut tubular elastic dressing retainer that is pre-shaped to accommodate the more difficult to-fit parts of the body. Ideal for securing hot or cold packs, dressings, and electrodes.

Product	Size	Applications	Pkg Unit/ Case
GL750	Stress Vest	Small/Medium	20/box
GL752	Stress Vest	Large/X-Large	20/box
GL600	Head	Universal Fit	20/box
GL620	Perineum Panty	Small/Medium	20/box
GL622	Perineum Panty	Large/X-Large	20/box
GL641	Knee, Foot, Elbow, Hand	Universal Fit	20/box

## Elastic Net Dressing (Retainer/Tube Gauze)



### Surgitube<sup>®</sup> Tubular Gauze (For use without applicator)

Tube gauze that secures dressings without tape. Allows maximum airflow and will not unravel when cut. Made of 100% Cotton. Does not contain natural rubber latex.

Product	Size	Description	Width	Pkg Unit/Case
<b>50 YARD ROLLS</b>				
GL241	1P	Small fingers, toes (WHITE)	5/8"	1/box
GL246	1P	Small fingers, toes (BEIGE)	5/8"	1/box
GL242	2P	Large fingers, toes (WHITE)	7/8"	1/box
GL247	2P	Large fingers, toes (BEIGE)	7/8"	1/box
GL242A	2AP	Infant legs, bulky fingers (WHITE)	1 1/8"	1/box
GL243	3P	Hands, feet, lower legs, shoulders (WHITE)	1 1/2"	1/box
GL248	3P	Hands, feet, lower legs, shoulders (BEIGE)	1 1/2"	1/box
GL244	4P	Wrists, elbows (tight weave) (WHITE)	1 1/2"	1/box
GL245	5P	Head, shoulders, breasts, legs (WHITE)	2"	1/box
GL260	5P	Head shoulders, breasts, legs (BEIGE)	2"	1/box
<b>15 YARD ROLLS</b>				
GL251	1	Small fingers, toes (WHITE)	5/8"	1/box
GL252	2	Large fingers, toes (WHITE)	7/8"	1/box
GL252A	2A	Infant legs, bulky fingers (WHITE)	1 1/8"	1/box
<b>10 YARD ROLLS</b>				
GL110W	1	Small fingers, toes (WHITE)	5/8"	1/box
GL210W	2	Large fingers, toes (WHITE)	7/8"	1/box
GL210F	2	Large fingers, toes (BEIGE)	7/8"	1/box
GL310W	3	Hands, feet, lower legs, shoulders (WHITE)	1 1/2"	1/box

Continued on next page

## Elastic Net Dressing (Retainer/Tube Gauze)



### Surgitube<sup>®</sup> Tubular Gauze (For use without applicator)

Product	Size	Description	Width	Pkg Unit/Case
<b>5 YARD ROLLS</b>				
GL105W	1	Small fingers, toes (WHITE)	5/8"	1/box
GL205W	2	Large fingers, toes (WHITE)	5/8"	1/box
GL205F	2	Large fingers, toes (BEIGE)	7/8"	1/box
GL305W	3	Hands, feet, lower legs, shoulders (WHITE)	7/8"	1/box
GL305F	3	Hands, feet, lower legs, shoulders (WHITE)	1 1/8"	1/box
GL405W	4	Wrists, elbows (tight weave) (WHITE)	1 1/2"	1/box



### SurgiGrip Tubular Elastic Support Bandage

Machine washable tubular elastic bandage designed to provide support and compression. Made of a soft cotton knit interlaced with elastic threads. Perfect for sprains and strains support following cast removal, soft-tissue injuries, joints and muscles, post-graft pressure for burn patients, and edema management/control. Can be applied to sensitive areas using an applicator and stays in one place with no clips or tape. 8-12 mm Hg compression at 50% stretch. 11 yard (10 meter) length. Does not contain natural rubber latex.

Product	Size	Applications	Pkg Unit/Case
GLB10	2 1/2" (6.25 cm)	Small hands and limbs	1/box
GLC10	2 3/4" (6.75 cm)	Adult hands, arms and legs	1/box
GLD10	3" (7.5 cm)	Large arms and legs	1/box
GLE10	3 1/2" (8.75 cm)	Legs or small thighs	1/box
GLF10	4" (10.00 cm)	Large knees or thighs	1/box
GLG10	4 1/2" (12.00 cm)	Large thighs	1/box
GLJ10	6 3/4" (17.5cm)	Small trunk	1/box
GLK10	8" (20.32 cm)	Medium trunk	1/box
GLL10	12 3/4" (32.38 cm)	Large trunk	1/box

## Elastic Net Dressing (Retainer/Tube Gauze)



### SurgiGrip Tubular Elastic Support Bandage

A tubular elastic bandage, containing latex, designed to provide support and compression for sprains and strains following cast removal, soft-tissue injuries, joints and muscles, post-graft pressure for burn patients, and edema management/control. Convenient alternative to flat wrap bandages. Can be applied to sensitive areas using an applicator and stays in one place with no clips or tape. Machine washable. 11 yard (10 meter) length.

Product	Size	Application	Pkg Unit/Case
WMB10	2 1/2" (6.25 cm)	Small hands and limbs	1/box
WMC10	2 3/4" (6.75 cm)	Adult hands, arms and legs	1/box



### Surgitube and SurgiGrip Tubular Gauze Applicators

For use with certain styles of Surgitube Tubular Gauze, so primary bandage is not dislodged. Also used with certain styles of SurgiGrip Tubular Elastic Support Bandages to eliminate possible irritation to sensitive areas



Product	Size	Applications	Pkg Unit/Case
<b>METAL CAGE APPLICATORS</b>			
GL230	0M	GL209, GL219	1/box
GL231	1M	GL209, GL219	1/box
GL232	2M	GL210, GL220, GLB10, WMB10, CL211, GL221, GLC10, GLD10, CLE10	1/box
GL234	3M	WMC10, WMD10, WME10, GL212, GL222, GLF10, GLG10, WMF10	1/box
GL236	4M	WMG10	1/box
<b>LUCITE TUBE APPLICATORS</b>			
GL231P	1L	GL209, GL219	1/box
GL232P	2L	GL210, GL220, GL211, GLC10, GLD10, GLF10	1/box
GL234P	3L	WMC10, WMD10, WME10	1/box
<b>ALUMINUM SPLINT TYPE APPLICATORS</b>			
GLALUAPP		GL209, CL219, GLn110W, GL210W, GL210F, GL105W, GL105F, GL205W, CL205F	12/case
<b>APPLICATORS SETS</b>			
GL227	1.2.3	Metal Cage Applicators	1 set/box

## Tapes/ Laparotomy Sponges



### Gentell FixTape

Gentell FixTape is a self-adhesive, non-woven fabric tape for the fixation of dressings, tubes and catheters with easy-to-remove protection paper. FixTape is latex free.

Product	Description	Pkg Unit/Case	HCPCS
GEN-10600	2" x 11yards	16/case	A4450
GEN-10640	4" x 11yards	16/case	A4450



### Ultrafix Self-Adhesive Dressing Retention Tape

A non-sterile and non-woven fixative tape that is breathable and especially gentle to the skin. Available in multiple widths that can be cut to size. Used to secure primary dressings such as gauze and foams. Does not contain natural rubber latex.

Product	Description	Pkg Unit/Case
62037	4" x 11yards	10/case
62038	6" x 11yards	10/case



### Laparotomy Sponges, Prewashed and X-Ray Detectable

These sponges are made of 100% cotton gauze, soft, absorbent, and are packaged pre-washed. All edges are folded in and stitched with cotton thread,

Product	Description	Pkg Unit/Case
1205	12" x 12", 4 ply, (molded pouch)	10/pkg, 40 pkgs/case, 400/case
1206	12" x 12", 4 ply, (flat pouch)	5/pkg, 60 pkgs/case,300/case
1213	18" x 18", 4 ply, (flat pouch)	5/pkg, 40 pkgs/case,200/case
1214	18" x 18", 4 ply, (molded pouch)	5/pkg, 40 pkgs/case,200/case
1223	4" x 18", 4 ply, (molded pouch)	5/pkg, 60 pkgs/case,300/case
12731	8" x 24 ", 8 ply, (molded pouch)	2/pkg, 50 pkgs/case,100/case
12751	8" x 24", 8 ply, (molded pouch)	1/pkg, 100 pkgs/case,100/case
12931	8" x 36", 4 ply, (molded pouch)	2/pkg, 50 pkgs/case, 100/case

## Compression Bandages & Bandaging Systems



### Sorban Undercast Padding Bandage

This bandage is a lightweight 100% Cotton bandage that offers protection of bony prominences and absorption of wound exudate. It is primarily used under other forms of compression bandaging and hard casts. Does not contain natural rubber latex.

Product	Description	Pkg Unit/Case
79504	4" x 5 yards	12/case



### Dufore Four-Layer Bandaging System

This system is designed to sustain effective levels of graduated compression for up to 7 days. The first layer, an absorptive padding layer, provides padding for bony prominences, distributes compression evenly over the leg and absorbs exudate. The second layer provides short-stretch compression by resisting expansion of the calf muscle during limb movement. It works with the third and fourth layers elastic compression bandages to effectively minimize edema and increase venous return. A wound contact dressing is included with the system to help maintain a moist wound healing environment and to protect the wound from outside contaminants. Does not contain natural rubber latex.

Product	Description	Pkg Unit/Case
72414	4 layer system	8/case



### Primer<sup>®</sup> Modified Unna Boot

A non-raveling thread locked gauze bandage that has been evenly impregnated with a non-hardening paste of zinc oxide, available with or without calamine. The Unna Boot is for use in a compression dressing for the lower leg. Ideal for managing extremities (e.g., varicose ulcers) requiring protection and support.

Product	Size	Applications	Pkg Unit/Case
GL3001	3" x 10 yards	12/case	1/box
GL4001	4" x 10 yards	12/case	1/box
GL3001C	3" x 10 yards (with Calamine)	12/case	1/box
GL4001C	4" x 10 yards (with Calamine)	12/case	1/box

## Compression Bandages & Bandaging Systems



### Unna-Pak: Primer<sup>®</sup> Modified Unna Boot and Duban Self Adherent Bandage

For the ambulatory management of venous stasis leg ulcers. A convenient system containing two bandages that together form a two-layer Unna Boot. The Unna-Pak Primer has all the advantages of the Primer<sup>®</sup> Modified Unna Boot with the additional benefit of a flesh colored secondary dressing that is lightweight and easy to apply.

Product	Size	Applications	Pkg Unit/Case
GL2003	3"	12/case	1/box
GL4001	4"	12/case	1/box



## Protective Products



### Glen-Sleeve II Arm and Leg Protectors

Knitted tubular sleeves designed to provide mild compression and protect the fragile skin of long-term care residents against everyday scrapes, shearing and bumps. Cotton and lycra construction with a polypropylene fiber that creates a wicking action to remove excess moisture from the skin. Durable sewn lockstitch construction with finished hems on both ends. Does not contain

Product	Description	Size	Pkg Unit/Case
GL1000	Hand, Wrist, Arm (WHITE)	18" (L) x 3" (W)	12 pairs/case
GL1000B	Hand, Wrist, Arm (BEIGE)	18" (L) x 3" (W)	12 pairs/case
GL1000WP	Padded Hand, Wrist, Arm (WHITE)	13 1/2" (L) x 3" (W)	12 pairs/case
GL1100B	Petite Hand, Wrist, Arm (BEIGE)	16 1/2" (L) x 2" (W)	12 pairs/case
GL2000	Hand, Wrist, Thumb, Arm (WHITE)	18" (L) x 3" (W)	12 pairs/case
GL3000	Below Knee (WHITE)	15" (L) x 3 1/2" (W)	12 pairs/case
GL3000B	Below Knee (BEIGE)	18" (L) x 3 1/2" (W)	12 pairs/case



### Glen-Sleeve Arm Protectors

Knitted tubular sleeves designed to provide mild compression and protect the fragile skin of long-term care residents against everyday scrapes, shearing and bumps. Made from polyester, cotton, nylon, and spandex.

Product	Description	Size	Pkg Unit/Case
GL899	Small/Medium (yellow band)	18" - 20" (L) x 3 1/4" (W)	12 pairs/case
GL900	Large/X-Large (blue band)	18" - 20" (L) x 3 1/2" (W)	12 pairs/case



### Heel & Elbow Protectors

Designed to protect against decubitus ulcers and other skin damage common among bed or wheelchair-bound patients. Consists of a knitted tube of soft fabric with a tricot-covered urethane foam pad stitched inside. Outer fabric moves independently of foam pad so the foam pad stays in contact with the

Product	Description	Size	Pkg Unit/Case
GL800	Small (yellow band)	9 1/2"	12 pairs/case
GL801	Medium/Large (blue band)	11"	12 pairs/case
GL802	X-Large (green band)	11 3/4" - 12 1/4"	12 pairs/case

## Gauze Burn Dressings



### Dupress Burn Dressings

These lint free pads help minimize trauma and adherence to the wound. The dressings are made of a fine mesh gauze cover with cellulose fill, to allow the dressing absorb fluid. Stitch construction with finished hems on both ends. Does not contain natural rubber latex.

Product	Description	Pkg Unit/Case
STERILE		
11903	4" x 4", 32 ply, Burn Gauze	10/pkg, 50 pkgs/case 20x12 mesh wrapped 500/case
11904	16" x 16", 16 ply, Burn Gauze	2/pkg, 40 pkgs/case, 20 x 12 mesh wrapped 80/case
11907	18" x 22" compress, rayon	1/pkg, 30 pkgs/case, cellulose wrapped in Gauze 30/case
11908	24" x 36" compress, rayon	1/pkg, 24 pkgs/case, cellulose wrapped in Gauze 24/case
8007	9" x 5 yds. Durlix, fluff gauze	1/pkg, 50 pkgs/case, 50/case
8136S	18" x 72" folded gauze, 20 x 12 mesh	2/pkg, 50 pkgs/case, 100/case
9825	9" x 9", 8 ply, gauze pad (10's )	10/pkg, 40 pkgs/case, 400/case

## Fixation Devices



### Cath-Strip<sup>®</sup> Reclosable Catheter Fastener

Cath-Strip<sup>®</sup> is a multi-use, reclosable catheter fastener. The material features flexibility and dynamic adherence which provide excellent tape-to-skin adhesion while minimizing the tension that may result in irritation, blistering, or skin shear. Cath-Strip is designed to accommodate a variety of catheter fastening requirements including; urinary catheters, gastrostomy and jejunostomy tubes, wound drainage systems, central line catheters, and multi-port IVs. Does not contain natural rubber latex stitch construction with finished hems on both ends. Does not contain natural rubber latex.

Product	Description	Pkg Unit/Case
CS50	1 pouch, 50/box	50/box



### UC Strip Catheter Tubing Fastener

UC Strip is an economical, one-piece catheter/tubing fastener. This flexible adhesive fastener secures urinary and gastrostomy catheter tubing to the patient without the irritation or tissue trauma caused by other tapes. UC Strip is made of a flexible material with a skin-friendly adhesive. The fastener allows limited side-to-side movement while it prevents the tubing from being pulled out. Does not contain natural rubber latex.

Product	Description	Pkg Unit/Case
UC20	1 pouch, 20/box	20/box



### NG Strip<sup>®</sup> Nasal Tube Fastener

NG Strip<sup>®</sup> is a securement device for any nasogastric or nasal feeding tube. NG Strip eliminates the time consuming ritual of securing nasal tubing with strips of irritating tape. NG Strip provides excellent tape-to-skin adhesion while minimizing the tension that may result in irritation, blistering, or skin shear. NG Strip also minimizes the chance of irritation and tissue necrosis that are common with conventional tapes. Does not contain natural rubber latex.

Product	Description	Pkg Unit/Case
NG25	Adult - 1/pouch, 25/box	25/box
NG40	Pediatric - 10/pouch, 4/box	10/pkg, 4 pkgs/box, 40/box
NG50	Adult- 10/pouch, 5/box	10/pkg, 5 pkgs/box. 50/box
NG60	Infant - 10/pouch, 6/box	10/pkg, 6 pkgs/box. 60/box

## Fixation Devices



### The Marsupial Pouch Wound Drainage Container Holder

The Marsupial pouch is an adjustable belt that fits comfortably around the waste. The pouch is designed to comfortably hold cumbersome wound drainage containers. It allows free movement of the arms and hands. Ideal for patients undergoing breast, chest, or abdominal surgery. The belt and pouch are made of terry cloth.

Product	Description	Pkg Unit/Case
MAR-01	1 per pouch	1/pkg

## Wound Closure Strips

### Suture Strip<sup>®</sup> Plus Flexible Wound Closure Strips - Sterile



Suture Strip<sup>®</sup> Plus is a wound closure strip that is made of macroporous non-woven polyamide with a pressure sensitive adhesive. The inherent elasticity of the product allows the strips to elongate and yield to the movement of skin, preventing the concentration of forces at the proximal ends and reducing the chances of traction blisters or skin shear. The strips are sterile, water resistant and have a skin-friendly adhesive. It is for use as a primary wound closure of surface tissue, traumatic wound closure, and for instances when aesthetic results are of great importance (such as in plastic surgery). Does not contain natural rubber latex.

Product	Description	Pkg Unit/Case
TP1100	1/8" x 3"	5/envelope, 50 envelopes/box
TP1101	1/4" x 3"	3/envelope, 50 envelopes/box
TP1102	1/4" x 4"	10/envelope, 50 envelopes/box
TP1103	1/2" x 4"	6/envelope, 50 envelopes/box
TP1104	1/4" x 1 1/2"	6/envelope, 50 envelopes/box
TP1105	1" x 5"	4/envelope, 25 envelopes/box
9825	9" x 9", 8 ply gauze pad (10's)	10/pkg, 40 pkgs/case, 400/case


### Suture Strip<sup>®</sup> Wound Closure Strips - Non-Sterile



Suture Strip<sup>®</sup> provides superior wound security, but with slightly less flexibility than Suture Strip Plus. Does not contain natural rubber latex

Product	Description	Pkg Unit/Case
TS3100	1/8" x 3"	5/envelope, 50 envelopes/box
TS3101	1/4" x 3"	3/envelope, 50 envelopes/box
TS3102	1/4" x 4"	10/envelope, 50 envelopes/box
TS3103	1/2" x 4"	6/envelope, 50 envelopes/box
TS3104	1/4" x 1 1/2"	6/envelope, 50 envelopes/box

## Wound Closure Strips



### Shur Strip Wound Closure Strips - Sterile

Shur Strip sterile skin closures are made of a porous, non-woven backing coated with a pressure-sensitive, hypoallergenic adhesive. These strips are not elastic, keeping wound edges securely together to maximize wound healing. Does not contain natural rubber latex.

Product	Description	Pkg Unit/Case
DKC81118	1/4" x 3"	3/envelope, 50 envelopes/box
DKC81120	1/4" x 4"	10/envelope, 50 envelopes/box
DKC81121	1/2" x 4"	6/envelope, 50 envelopes/box
DKC81119	1/4" x 1 1/2"	6/envelope, 50 envelopes/box



### Shur Strip Wound Closure Strips Non-Sterile Bulk

These wound closure strips are bulk packaged and are a non-sterile version of our Shur Strip. Does not contain natural rubber latex

Product	Description	Pkg Unit/Case
DKC71153	1/2" x 4"	6/card, 11 cards, 6600 total





Gentell

2701 Bartram Road, Bristol, PA

19007

[www.gentell.com](http://www.gentell.com)

(215) 788-2700

Gentell is one of the largest vertically-integrated wound care companies in the world. We provide efficient, affordable, patient-specific wound treatments to patients in nursing homes, hospices, hospitals and other settings. We employ clinical specialists throughout the United States and Canada to assist with treatment recommendations and provide education about the proper utilization of wound care products. With production facilities in the United States, Canada, Mexico and China, we are one of North America's largest manufacturers of wound care supplies.